AIR COMPRESSION THERAPY SYSTEM

MODELLO 4C DIGITAL/APT 5004

PRESSOTERAPIA
ITALIA

DINAMICAMED

Reg. Ministero della Salute

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Dear Customer,

We thank you for choosing Air Compression Therapy System,
designed with high quality materials.

For correct use, we remind you to read the manual before use.

Index.....	p. 22
Pressure Therapy and its applications.....	p. 23
Components.....	p. 25
Control panel	p. 26
Before use	p. 27
Massage Programmes	p.30
After use	p.31
Cleaning and maintenance	p.32
Warnings and suggestions	p.34
Contraindications.....	p.38
Warranty	p.39
Technical features	p.41

Pressure Therapy is a treatment that improves the functioning of the circulatory system and the lymphatic system. It consists of applying external pressure to the area to be treated, using a machine specially designed for this purpose. It is a simulation of manual lymphatic drainage that improves the functioning of the circulatory system and the lymphatic system and helps the natural processes of elimination from the body of toxic substances.

The lymphatic drainage has the purpose of reducing water retention and promoting the flow of lymph inside its vessels, thus avoiding stagnation in some areas of the body.

The pressure therapy is indicated for:

- Reduce the unsightly appearance of cellulite.
- Reduce water retention, caused by the stagnation of lymph in one or more sites in the human body (both congenital and post-operative) which causes pain and swelling.

- Tone the muscles
- Improve lymphatic and venous circulation.
- Help the elimination of toxic substances from the body.
- Promote the slimming of some parts of the body, for example legs.
- Reduce traumatic edema.
- Reduce abdominal bloating.
- Relieve muscle fatigue and discomfort caused by swollen, heavy or tired feet, ankles and legs.

Usually, the anatomical areas most treated with the pressure therapy are lower limbs, upper limbs, abdomen and buttocks. Generally, leg treatments are the most popular. The air compression therapy system is composed of an electronic device connected to cushions that are applied around the anatomical area to be treated.

Through a system of air jets, the machine exerts an intermittent and sequential pressure on the treated area, favoring the movement of both lymphatic and venous flows following the natural path that goes from the periphery of the limb to the heart.


OPTIONAL ACCESSORIES

ARM CUFF


WAIST CUFF



NOTE: When you receive the machine, first time before starting to use the unit, please make sure open the filter cover, and remove the security screw. The screw **ONLY** are used for fixing and protecting the control unit from damage during transportation.


STEP 1 Insert the main plug of power line into electrical outlet of 110V~120V or 220~240V

STEP 2 Check if the hose has been connected well to the cuff(s): if yes, skip to STEP 3; if not, please fit the pins of the hose into the outlets of cuff(s).


STEP 3 Connect the main end of the connector hose to the machine

STEP 4 Check and ensure all connections to the Control Unit and cuffs make a good seal.

STEP 5 Put on the cuffs

NOTE: The machine had better be powered up prior to attaching cuff(s)

STEP 6 Press the power button. 

STEP 7 Set the pressure by pushing the +/- button  

Note: for first time use, it had better start with the pressure not too high, then could adjust the pressure to the level that makes you feel comfortable.

STEP 8 Set the time by pushing the  button

Note: time range is 1~60 minutes. Press - / + to subtract/ add 10 minutes.

STEP 9 Choose the mode by pushing button

STEP 10 Set the position by pushing


Note: The original setting is light on, means all chambers will work; If you want to skip one position like 1, press it close the chamber and light off, means the chamber won't work. If need restart the chamber, press again to light on. When inflating or deflating each chamber, that chamber light is flashing.


STEP 11 Press the start/ pause  button

Note: Original state is pause, the light is green; Press it to start, the light turns into blue.

*** DO NOT unzip the cuff(s) when the unit is switched on. This may cause damage to the cuff(s). Always turn off the power before unzipping the garment.**


A). Lymphatic drainage: contrasts orange peel skin. The pressure is exerted by only one chamber at a time, starting from the foot up to the groin.


B). Deep lymphatic drainage: toning and anti-cellulite. The pressure is exerted evenly and continuously. The chambers are all inflated progressively


C). Beauty: massage 2+2

D). Total Relax: total relaxation session. The pressure comes exercised on all chambers simultaneously.


E). Slimming-Toning: acts by toning the legs.


F). Anti-cellulite: program that contrasts orange peel skin + toning.

- Just needs to disconnect the hose from the machine by pulling the main end out


- The cuffs connected with hose are always kept and used as a whole; When use, it just need insert the other end of the hose the machine


- Had better not disconnect the hose from the cuffs frequently; Abuse may cause damage to the pins or hose.


Note: When travel with the product, it could disconnect the hose from the cuffs to save space.

No pulling the hose


Disconnect the hose from cuff by holding pins to pull up.

Do not use water, oil, benzene, alcohol, gasoline and chemical agents, to wash the machine.

· Otherwise, the lifetime of the machine will be shortened.

Please wipe it with a dry cloth.


Do not use water, oil, benzene, alcohol, gasoline and chemical agents, to wash the garments.

· Otherwise, the lifetime of the garments will be shortened.

Please wipe garments with damp cloth.


Do not wash your cuffs, rub with a towel gently.

NOTE: All cuffs are for air pressure massage, uneasy to repair.

Please be careful to use and keep!

Do not use water, oil, benzene, alcohol, gasoline and chemical agents, to wash the machine.


Otherwise, the lifetime of the machine will be shortened.

Please wipe it with a dry cloth.

For long-time keeping, please put the equipment in the box.

The environment for use and store

At room temperature, in a clean and dry place to use and store the equipment. (Temperature: 10°C to 40 °C, humidity: 30% to 85%).

- Please read the guide carefully before using the product
- Please place the guide closed to the machine to read it better
- In order to avoid danger or some non-standard operations, to ensure safety, here the guide lists some details, and the users had better comply with them when using the product.

If feel abnormal or painful when using,
please stop using the product immediately.


Do not use the machine for nudity.
Wear underwear/ clothes before wear the garment.

For first time, it had better not exceed 20 minutes.

For first time to use the machine, it had better last about 10 minutes, and then according to user's own physical situation, and gradually extended treating time.

When using the machine, do not go with
the cuffs.

Do not insert a few power into the same socket at the same time, it had better plug ONLY a power to a socket.

- May cause fire, electric shock, or other damage.

Keep the hand dry when plugging in or unplugging the power supply.

- May cause leakage, electric shock

Only use 220V~240V or 110V~120V single-phase current, after using, please unplug the power correctly. • May cause fire, electric shock, or other damage.


Do not bend and fold the power cord.


Power line can not bear the weight, can not be tied in the tables and chairs.

- May cause leakage, electric shock, or fire.

Do not hit or beat the machine.


- Such abuse will damage the machine.

Do not put any foreign item on the equipment.

- If the sensor is clogged, placing the item on the conduit or bending the power cord, it will cause a pressure disorder.

(the equipment may be damaged or abnormal)

Be careful! Do not put nails, needles, staplers, metal or food into the machine.

- May cause short circuit or conduction.

Use the equipment for massage only.

- Abusing the machine may cause damage.


IT IS NOT A TOY FOR KIDS AND PETS.

Prohibit being used by infants and young children

· May cause damage such as electric shock


The equipment can not be repaired and maintained during use. Be sure to perform maintenance when power off.

Protect the equipment from direct sunlight or heat, such as gas stove.

The equipment will be deformed or discolored.

When the equipment, button cell, cuffs and package carton box are discarded, it should be disposed of in accordance with the garbage disposal procedures.

Pay attention to environmental pollution.

Do not use the machine even if it is alone one of the following diseases:

- General ill feeling, flu or illness
- Renal / cardiac / liver insufficiency
- Arterial insufficiency
- People with pacemakers
- Injury /Acute dermal inflammation
- Dermatological disorders (e.g. lymphangitis)
- Infections
- Gangrene or dermatitis
- Arrhythmia
- Plexopathy and neuropathy
- Severe arteriosclerosis
- Ischemic vascular diseases
- Deformity of the limbs or presence of metal joint prostheses
- Deep thrombophlebitis or varices
- Pulmonary edema
- Pulmonary embolism
- Wounds or injuries
- Cancers in the area of application or the near area
- Unstable hypertension
- Diabetes mellitus
- Cirrhosis of the liver

It is forbidden to use the machine even for pregnant women.

ATTENTION! Do not use when in doubt, except before you have consulted your doctor.

- The device is guaranteed to the first user for the period of twenty-four (24) months, three (3) months for applicators and accessories from the date of purchase against defects in materials or workmanship, as long as used appropriately and maintained in normal conditions of efficiency.
- The validity of the guarantee, in the event that the user uses it for professional purposes, is twelve (12) months from the date of purchase against material or manufacturing defects, provided it is used appropriately and maintained in normal conditions of efficiency.

The validity of the guarantee is limited in the following cases:

- No guarantee for accessories and materials to be considered "consumables".
- This Warranty is valid and effective in the country where the Product has been purchased.
- In the event that the product was purchased in any country of the European Community, the guarantee is still valid in all its countries.
- Problems due to drops or incorrect use of voltage.
- If the device has been opened or modified by unauthorized personnel.
- If the appliance has been used in unsuitable environments.
- If the appliance has not been used following the instructions of this manual

Equipment warranty: 2 years;

Warranty on accessories and applicators: 3 months.

To use the warranty service, the user must comply with the following warranty clauses:

- The products must be delivered for repairs by and at the expense of the Customer in their original packaging.
- The product warranty will be subject to the presentation of a fiscal document (receipt or sales invoice), which certifies the date of purchase of the product.
- The repair will not affect the original expiration date of the guarantee and will not lead to its renewal or termination.
- If, at the time of repair, no defect is found in the product itself, the costs relating to the technical verification will still be charged.
- The warranty is void if the failure has been caused by bumps, falls, incorrect or improper use of the product, use of a power supply or non-original charger, accidental events, product alteration, replacement / detachment of warranty seals and / or product tampering.
- Furthermore, it does not cover damage caused during transport due to the use of unsuitable packaging.
- The warranty is not liable for the impossibility of using the product, for other incidental or consequential costs or for other expenses incurred by the buyer. Furthermore, it does not cover damage caused during transport due to the use of unsuitable packaging.
- The warranty is not liable for the impossibility of using the product, for other incidental or consequential costs or for other expenses incurred by the buyer.

CUSTOMER CARE

In the event of a breakdown or malfunction of the device, contact your dealer.

THE PACKAGING CONTAINS

- Machine body with power supply
- Rubber hose with connector
- Cuffs
- User manual

TECHNICAL FEATURES

Trade name: Air Compression Therapy System

Model: 4C DIGITAL/APT 5004

Adjustable treatment time: from 1 to 60 minutes

Adjustable pressure: from 30 to 250 mmHg

Power supply: 220-230V 50 / 60Hz 65VA

To receive assistance or for any further information, contact our Customer Service:

+39 3939275259

+39 3939337057

Made in PRC

CLASS IIa MEDICAL DEVICE

Imported by: NITOR MEDIA SRL

Via Montenapoleone 8, 20121 MILANO (MI),

in compliance with Directive 93/42 / EEC

CE Certificate N. 0598 / FI21 / 07003.

The Serial Number label is positioned under the device

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Handelsname: Air Compression Therapy System

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the 1990s, the number of people in the world who are under 15 years of age has increased by 1.2 billion, from 1.1 billion in 1980 to 2.3 billion in 1999. The number of children under 15 years of age in the world is projected to increase to 3.1 billion by 2015, with the largest increases occurring in the developing countries (United Nations, 1999).

There is a growing awareness of the need to address the needs of children in the world, and the United Nations has established the Convention on the Rights of the Child (CRC) to protect and promote the rights of children. The CRC is a legally binding treaty that sets out the rights of children and the obligations of governments to protect and promote these rights. The CRC has been ratified by 112 countries, and it is the most widely ratified human rights treaty in the world (United Nations, 1999).

The CRC is a comprehensive document that covers a wide range of issues, including the right to life, the right to health, the right to education, the right to play, and the right to be heard. The CRC also sets out the obligations of governments to protect and promote the rights of children, and it provides a framework for monitoring and evaluating the implementation of the CRC.

The CRC is a landmark document in the history of children's rights, and it has played a major role in the development of children's rights around the world. The CRC has inspired the development of national laws and policies that protect and promote the rights of children, and it has provided a framework for monitoring and evaluating the implementation of these laws and policies.

The CRC is a comprehensive document that covers a wide range of issues, including the right to life, the right to health, the right to education, the right to play, and the right to be heard. The CRC also sets out the obligations of governments to protect and promote the rights of children, and it provides a framework for monitoring and evaluating the implementation of the CRC.

The CRC is a landmark document in the history of children's rights, and it has played a major role in the development of children's rights around the world. The CRC has inspired the development of national laws and policies that protect and promote the rights of children, and it has provided a framework for monitoring and evaluating the implementation of these laws and policies.

The CRC is a comprehensive document that covers a wide range of issues, including the right to life, the right to health, the right to education, the right to play, and the right to be heard. The CRC also sets out the obligations of governments to protect and promote the rights of children, and it provides a framework for monitoring and evaluating the implementation of the CRC.

The CRC is a landmark document in the history of children's rights, and it has played a major role in the development of children's rights around the world. The CRC has inspired the development of national laws and policies that protect and promote the rights of children, and it has provided a framework for monitoring and evaluating the implementation of these laws and policies.

The CRC is a comprehensive document that covers a wide range of issues, including the right to life, the right to health, the right to education, the right to play, and the right to be heard. The CRC also sets out the obligations of governments to protect and promote the rights of children, and it provides a framework for monitoring and evaluating the implementation of the CRC.